# PRB 1026 BIOILOGY III Laboratory Manual Semester 2 SECOND





Arkih QH 317 P921 2016

Edited by: Mohd Ridwan Abd Rahman Muhamad Ikhwan Idris Rohaiza Daud



UNIVERSITE MALAYSIA SARAWAK

•

## PRB 1026

.

## BIOLOGGY JU Laboratory Manual Semester 2



•

#### Pusat Khidmat Maklumai Akademik UNIVERSPERMALAYSI SARAWAK

**\_\_\_** 

P.KHIDMAT MAKLUMAT AKADEMIK UNIMAS



## PRB 1026 BIOLOGGY II Laboratory Manual Semester 2

---



-

Mohd Ridwan Abd Rahman

Muhamad Ikhwan Idris

Rohaiza Daud

Mohd Razif Othman

Roberta Chaya Tawie Tingga

Christharina S Gintoron

Mohamad Aminudin Mustapha

Mohamad Fhaizal Mohamad Bukhori

Universiti Malaysia Sarawak

Second Edition 2016
Third Publication 2015
Second Publication 2014
First Publication 2013
© Mohd Ridwan Abd Rahman,
Muhamad Ikhwan Idris.

Rohaiza Daud, Mohamad Razif Othman, Roberta Chaya Tawie Tingga, Christharina S Gintoron, Mohd Aminudin Mustapha, Mohamad Fhaizal Mohamad Bukhori,

All rights reserved. No part of this publication may be reproduced, stored in retrieval system, or transmitted, in any form or by any means, electronic, mechanical, photocopying, recording or otherwise, without the prior permission of the publisher.

.

Published in Malaysia by UNIMAS Publisher,

Universiti Malaysia Sarawak, 94300 Kota Samarahan,

Sarawak, Malaysia.

Printed in Malaysia by Malien Press Sendirian Berhad G/Floor Lot 184, Section 49, Jalan Abell, 93100 Kuching, Sarawak, Malaysia.

Perpustakaan Negara Malaysia

Cataloguing in-Publication Data

PRB 1026 BIOLOGY ii : Laboratory Manual. Semester 2 / Editors Name Mohd Ridwan Abd. Rahman, Muhamad Ikhwan Idris, Rohaiza Daud, Roberta Chaya Tawie Tingga, Mohamad Razif Othman, Christharina S Gintoron, Mohd Aminudin Mustapha, Mohamad Fhaizal Mohamad Bukhori, ; Contributors Name Ramlah Zainuddin, Siti Akmar Khadijah Ab Rahim Bibliography: pages 71-72 ISBN 978-967-2008-19-4
1. Biology--Laboratory manuals.
1. Mohd. Ridwan Abd Rahman, II. Muhamad Ikhwan Idris.
IH. Rohaiza Daud, Iv. Mohamad Razif Othman.
V. Roberta Chaya Tawie Tingga, VI. Christharina S Gintoron.
VII Mohd Aminudin Mustapha, VII. Mohamad Fhaizal Bukhori, 1978-.
570.78

Artib QH 317 2016 Pusat Khidmat Maklumat Akademik UNIVERSITI MALAYSIA SARAWAK

## TABLE OF CONTENTS



vii

.

PREFACE	ix
ABBREVIATIONS	xi-xii
LABORATORY REGULATIONS	xiii-xiv
LABORATORY INSTRUCTIONS	xv-xvi
<b>Proctical 1</b> . DNA Extraction using Cetyltrimetyl Ammonium Bromide	1-7

: DNA Extraction using Cetytermetyr Annonnun Dronnue Practical 1 **T** 1 (CTAB)

: DNA Amplification using Polymerase Chain Reaction (PCR) 9-16 **Practical 2** 

Practical 3	•	Animal Evolution	17-28
Practical 4	•	Plant Taxonomy and Classification	29-42
Practical 5	•	Animal Identification	43-52
Practical 6	•	Water Quality Observation	53-64
Video Presentation and Reflection – "The Walking Whale"		65-69	
REFERENCE	S		71-72

٤

### ACKNOWLEDGEMENT

Special thanks and appreciation to all members of the Biology Unit in *Pusat Pengajian Pra-Universiti* (PPPU), Universiti Malaysia Sarawak (UNIMAS); Mohd Ridwan Abd Rahman, Muhamad Ikhwan Idris, Christharina S Gintoron, Mohamad Razif Othman, Mohd Aminudin Mustapha, Roberta Chaya Tawie Tingga, Rohaiza Daud and Mohamad Fhaizal Mohamad Bukhori who had been tirelessly involved in initiating, editing and compiling this laboratory manuals in order to make the second publication of this

manual a reality.

We also wish to thank the following people who helped in contributing the ideas for the second edition of this laboratory manual:

AP Dr Ramlah Zainudin, Faculty of Resource Science and Technology; AP Dr Muna Sabri, Faculty of Medical and Health Sciences; Prof. Dr. Lela Hj. Su'ut, Faculty of Medical and Health Sciences; Dr Siti Akmar Khadijah Ab. Rahim, Faculty of Resource Science and Technology; and Prof. Dr. Ahmad Hata Rasit, Faculty of Medical and Health Sciences.

Finally, we also like to thank Mr. Mohamad Nazri Mario from Publication Division, UNIMAS for contributing the book design and photo of the cover page for this PRB1026 Biology II Lab Manual Semester 2 (Second Edition).

## PREFACE

This second edition of *Biology Laboratory Manual Semester II* is published to meet the needs of the students taking Biology course in Pusat Pengajian Pra Universiti (PPPU), Universiti Malaysia Sarawak (UNIMAS). The purpose of this manual is to provide the students with an organised, user-friendly manual to aid their understanding on the practical aspects of biology as well as enhancing learning and understanding of the lectures. This manual has been published based on reviews from PPPU Biology lecturers and advice from the experts consisted of experienced lecturers teaching in degree courses in UNIMAS.

Each section of the manual has been designed according to the standard of a scientific report format. In **Introduction**, brief explanation is provided to help the students form better idea on the background of the laboratory before they can start with the practical. This is followed by **Objectives**, which states the aim of the laboratory and what should be achieved at the end of the laboratory session. In **Materials**, a list of materials including the instruments and apparatus used for each laboratory is prepared for the student's reference. In **Methodology**, step-by-step procedure is explained in detail to guide the students in carrying out the practical correctly. The **Results** section is where the students record important observations and findings that they have obtained. **Discussions** section is for the students to explain the results that they have obtained and to relate with their learning units, before summarizing it in the **Conclusions** 

section.

This manual is also loaded with a number of additional features. This includes the tearoff section to enable the students to submit their report at the end of laboratory session. Besides that, the topics of the manual have been arranged according to their weekly learning units in their syllabus. Notes and theories/operating principles are also included as students' guidance. Relevant diagrams are inserted to enhance the understanding of the description in the manual. Post-Lab Questions are also prepared not only to validate the student's understanding but also to enhance their understanding of the laboratory procedure and recognise its applications.

Therefore it is hoped that the publication of this *Biology Laboratory Manual Semester* II will serve as a useful guide for the students in performing their laboratory sessions.

## ABBREVIATIONS

%

-

percentage

°C	celsius
μl	microlitre
μS/cm	Microsiemens per centimetre
µmhos/cm	Micromhos per centimetre
AN	Ammoniacal nitrogen
BOD	Biochemical oxygen demand
CIA	Chloroform-isoamyl alcohol
cm	Centimetre
COD	Chemical oxygen demand
CTAB	Cetyltrimetylammonium bromide
Cyt c	Cytochrome c
ddH2O	Deionised distilled water
DNA	Deoxyribonucleic acid
dNTP	Deoxyribonucleotide triphosphate
DO	Dissolved oxygen
DOE	Department of Environment
eg	'exempli gratia' / for example
etc.	Et cetera
FNU	Fromazin nephelometric unit
FTU	Formazin turbidity Unit

GPS	<b>Global Positioning System</b>
i.e	'id est' / that is
KCl	Potassium chloride
mg/L@mg/l	Milligram per litre
$\mathbf{MgCl}_2$	Magnesium chloride
ml	Mililitre
mM	Milimolar

PRB 1026, Biology II Laboratory Manual Semester 2

.

NaCl	Sodium chloride
NH <sub>3</sub> N	Ammoniacal nitrogen
NTU	Nephelometric turbidity units
NWQS	National Water Quality Standards
PCR	Polymerase Chain Reaction
pН	Potential hydrogen
PPM	Parts-per-thousands
PPT	Parts-per-million
PSU	Practical salinity unit
rpm	revolutions per minute
sp.	Species
SS	Suspended solid
Taq	Thermus aquaticus
Tris – HCL	Tris(hydroxymethyl)aminomethane hydrochloride
UV	Ultraviolet
WQI	Water Quality Index



PRB 1026, Biology II Laboratory Manual Semester 2

.

### LABORATORY REGULATIONS

#### STUDENTS ARE STRICTLTY PROHIBITED TO ENTER LABORATORY

#### WITHOUT PERMISSION OF THE AUTHORIZED PERSONNELS (LECTURERS, TUTORS AND LAB DEMONSTRATORS)

- 1. Attendance to the laboratory is compulsory. A Medical Certificate is needed for
  - students who are unable to attend practical classes due to health reasons.
- 2. Students must be punctual to their practical sessions. Any latecomers will not be permitted to enter the laboratory.
- 3. Students must wear laboratory coat and proper/suitable shoes at all times. Sandals or open-toed shoes are strictly prohibited in the laboratory.
- 4. Read and follow all directions or methods of the experiment. Students should be prepared and understood all the procedures before performing the practical. Ask the lab demonstrators if there are any doubts about any part of the experiment.
- 5. Students are not allowed to eat and drink in the laboratory.
- 6. Students must wear safety goggles when working with chemicals, burners, or any substance that might get into the eyes.
- 7. Students must wear gloves whenever handling dangerous and corrosive chemicals.
- 8. Keep laboratory area clean and tidy. All instruments must be kept neat and clean at all times. DO NOT handle any equipment unless it is stated in the laboratory procedure.
- 9. Be serious and alert when working in the laboratory. Never mix chemicals just for fun because it might produce dangerous and explosive substances.

10. Do not waste any chemicals and take extreme care not to spill any material in the laboratory. Ask lab demonstrator about proper clean up procedure. Do not simply pour chemicals into the sink or trash container.

11. All accidents in the laboratory have to be reported to the lab demonstrators.

12. Take extreme care when handling scalpels or razor blades. Never point anything sharps towards other students. Students must immediately notify lab demonstrator if any student accidentally cut themselves or receive a cut.

13. Treat all living things with care and respect. Do not simply touch any organisms and specimens in the laboratory unless students have been given permission to do so. There are plants which are poisonous or with thorns and animals that may bite or scratched if alarmed.

14. Animals should be handled only if necessary. Lab demonstrator will give instructions on how to handle each species that may be brought into the laboratory.

15. Rinse any acids off the affected skin or clothing with water. Immediately inform lab demonstrator of any spillage.

16. Be aware of the location and proper use of safety equipment such as first aid-kit, fire extinguisher and emergency shower.

17. Students must clean and tidy the benches and return all equipment to its proper place when the experiment is completed. Wash hands after each experiment.

## STUDENTS WHO FAIL TO ABIDE BY THE ABOVE STATED RULES AND REGULATIONS ARE SUBJECTED TO STRICT DISCIPLINARY ACTION

Biology Coordinator, Pusat Pengajian Pra-Universiti, Universiti Malaysia Sarawak.



### LABORATORY INSTRUCTIONS

1. Students are advised to read the laboratory manual and make notes before attending each practical class in order to be adequately prepared. This will save the time and would ensure the students will be able to comprehend the practical in time.

- 2. Students are advised to bring the following items to every practical classes:
  - i. Laboratory Manual
  - ii. Any essential stationeries iii. Blank and Ruled A4 papers
- 3. Read and understand the practical guide in the laboratory manual before undertaking any practical exercise or experiment. Consult with the lecturers, demonstrators or lab assistant for any inquiry.
- 4. Submit a lab report to the lecturers/demonstrators at the end of each practical class in order to be marked.
- 5. The lab report depends on the nature of the practical: either experimental observations or drawings in plant or animal diversity topics.
- 6. Answer all questions and draw the observed specimens/slides/diagrams in the spaces provided.
- 7. Drawings or diagrams should be large (as a minimum of half a page), accurate, clear, titled and labeled accordingly. Scientific names of organisms should be written properly.
- 8. For drawings of specimens/slides observed under a microscope, the size of the objective must be given (X10, X40, etc.). Draw only what you visualized. YOU MUST NOT copy any drawings directly from other sources such as books or photocopied materials but students may use them as a guide.

#### SCIENTIFIC REPORT

\* If the practical class involves experiments, students have to write a complete scientific report in the following format;:

1. Title

2. Introduction

3. Materials and Methods

4. Results

5. Discussion

## Report should be in full sentences, past tense and comprehensible with precision, clarity and concise.

#### Title

The title must be short and precise.

#### Introduction

The introduction will be explaining the objectives of the experiment (maximum 100

words).

#### **Materials and Methods**

Summarize the procedures and the materials that were used.

#### Results

The data must be presented in tables, graphs or histograms whichever related, and with brief description.

#### Discussion

In the discussion, the results should be interpreted with valid explanation. In the case of obtained results are not as expected, students should discuss both the obtained and the expected results. Students may also make comparison on different methods and discuss on the related difficulties encountered.

#### Conclusion

The conclusion must be short and precise, based on the results and discussion

PRB 1026, Biology II Laboratory Manual Semester 2

•

**Practical 1 DNA Extraction Using Cetyltrimetylammonium bromide (CTAB)** 

Roberta Chaya Tawie Tingga and Mohd Ridwan Abd Rahman

#### Introduction

Deoxyribonucleic acid (DNA) isolation is an extraction process of DNA present in the nucleus of the cell. Methods used to isolate DNA are dependent on the source, age, and size of the sample. Isolation of DNA is needed for genetic analysis, which is used for scientific, medical, or forensic purposes. Common sources for DNA isolation include blood, hair, sperm, bones, nails, tissues, blood stains, saliva, epithelial cells, urine, bacteria, animal tissues, or plants. The isolation of DNA usually begins with lysis, or breakdown of tissue or cells. This process is essential for the destruction of protein structures and allows for release of nucleic acids from the nucleus. Lysis is carried out in a salt solution, containing detergents to denature proteins or proteases. In this experiment, the salt solution used is Cetyltrimetylammonium bromide (CTAB).

#### **Objectives**

- 1. To understand the process and principle of standard DNA isolation.
- 2. To extract DNA of animal tissue using Cetyltrimetylammonium bromide (CTAB) protocol.

#### Materials

#### 1. Animal tissue

- 2. Micropipettes
- 3. Microcentrifuge tubes
- 4. Dissecting kit
- 5. Aluminum foil
- 6. Water bath
- 7. Centrifuge machine
- 8. CTAB solution

Practical 1: DNA Extraction Using Cetyltrimetyl Ammonium Bromide (CTAB)

9. Proteinase K

10. Chloro-isoamyl alcohol (CIA)

11. Absolute ethanol (99.9% EtOH)

12.70% ethanol (70% EtOH)

13. Sodium chloride (3M NaCl)

14. Deionized distilled water (ddh<sub>9</sub>0)

#### Methods

- 1. Cut up and mince a small amount of animal tissue on the aluminum foil.
- 2. Then put 700  $\mu$ l of CTAB into 1.5 ml microcentrifuge tube.
- 3. Add 5  $\mu$ l of Proteinase K and incubate in water bath at 65°C for 20-30 minutes or longer until the tissue is completely dissolved.
- 4. After tissue has dissolved, add 700  $\mu$ l of CIA and inverse for 1 minute. Centrifuge at 13000 rpm for 10 minutes.
- 5. Take 550 µl of the upper layer of the supernatant and transfer to a newly labeled tube (do not touch the middle layer).
- 6. Add 550  $\mu$ l of cold absolute EtOH. Then, spin the tube at 13000 rpm for 10 minutes.
- 7. Discard EtOH by pouring off. Make sure that the pellet is still intact at the bottom of the tube. Observe if the pellet is visible (yellow or white coloured pellet).
- 8. Add 550  $\mu$ l of cold 70% EtOH, and 25  $\mu$ l of 3M NaCl. Spin at 13000 rpm for 10 minutes.
- 9. Discard EtOH. Observe if DNA pellet is still in the tube.
- 10. Re-dissolve the pellet in 25-30  $\mu$ l of (ddH<sub>a</sub>O).
- 11. Run 1-2  $\mu$ l of extraction product through 1% agarose gel electrophoresis. 12. Visualized the gel under UV transilluminator.

Name:





#### Laboratory Report

Write **TWO (2)** pages of handwritten report based on the DNA extraction experiment that you have conducted today. By using your own words, this report must include a brief **Introduction** (2 marks), **Objectives** (2 marks), **Methodology** (5 marks), Results (6 marks – including gel documentation), **Discussion** (10 marks), **Conclusion** (2 marks) and list of **References** used (3 marks). [30 marks]

Practical 1: DNA Extraction Using Cetyltrimetyl Ammonium Bromide (CTAB)



Name:	Matric No:	Group:	
Gel documentation	Pusat Khidmat Maktumat Akademik UNIVERSITI MALAYSIA SARAWAK	[3 marks]	
	·		

.

#### References

#### [3 marks]

Practical 1: DNA Extraction Using Cetyltrimetyl Ammonium Bromide (CTAB)

• •



٠



[3 marks]

#### **Post Laboratory Questions**

1. What is the role of CTAB solution in DNA isolation?

2. What does the Proteinase K do?

#### [1 mark]

3. What does the CIA do and why?

[4 marks]

-

4. How can we shorten the incubation time during the lysis process? [3 marks]

#### 5. What is the use of every centrifugal step in the protocol?

[1 mark]

6. Why does the DNA clumps together?

-

[4 marks]

•

.

Practical 1: DNA Extraction Using Cetyltrimetyl Ammonium Bromide (CTAB)

.

#### **Practical 2 DNA Amplification Using Polymerase** Chain Reaction (PCR)

Mohd Ridwan Abd Rahman and Roberta Chaya Tawie Tingga

#### Introduction

The polymerase chain reaction (PCR) is a technique to amplify the small amount of DNA up to a millions of copies by using the thermocycler. It involves enzymatic synthesis of the specific DNA including a combination between a DNA sample with the oligonucleotide primers, deoxynucleotide triphosphate (dNTPs) and the DNA polymerase in the buffer. There are three events that must take place in the PCR process; 1) template denaturation which is the most important process to make sure the cells' membrane are being lyses; 2) annealing of the primers and lastly; 3) the extension of the annealed primers using DNA Taq polymerase. The PCR technique is widely used in molecular biology, microbiology, genetics, diagnostics, clinical laboratories, forensic science, environmental science,

hereditary studies, paternity testing, and many other applications.

#### **Objectives**

- 1. To explain the function of each component used in the PCR amplification.
- 2. To explain the process involved in the amplification of DNA.
- 3. To give hands-on experience to the students on DNA amplification using the thermocycler.

#### NOTES

- Precautions should be taken to guard against contamination of the reaction with • trace amounts of DNA that could serve as templates.
- Disposable gloves should be worn throughout the preparation of the PCR. •
- Always include a control that contains all the components of the PCR except the template DNA.

Practical 2: DNA Amplification Using Polymerase Chain Reaction (PCR)

Matric No:

Group:

### Materials

- 1. 10x reaction buffer
- 2.  $MgCl_2$  (25 mM)
- 3. dNTPs mix (10mM)
- 4. Primers : \_\_\_\_\_\_ gene;
  - a. Forward \_
  - b. Reverse -
- 5. Template DNA
- 6. Taq DNA polymerase (5 u/µl)
- 7. Deionized distilled water (ddH<sub>0</sub>)
- 8. Thermocycler

#### Methods

1. First, make a master mix in the 1.5 ml microcentrifuge tube. Add each components in the following order:

#### Table 2.1: Components of PCR master mix.

No.	Component	V	×x	
		1x Reaction (25 µl)	y reactions	
1	10x reaction buffer	2.5		
2	$MgCl_2$ (25 mM)	1.5		
3	dNTPs mix (10 mM)	0.5	fos inov erus es	
4	Forward (10 mM)	1.0	sonkting sten	
5	Reverse (10 mM)	1.0		
6	ddH <sub>2</sub> O	**Z		
	Template DNA	a		
	Taq polymerase, (5 u/µl)	0.2		
	Total	25		

 $*_{x} = (v) (y)$ 

## \*\* $z = 25 \mu l \cdot \Sigma$ reaction components

- a = vary, depending on the quality and quantity of DNA Template
- $b = \sum$  Component reaction 1 until 6

Practical 2: DNA Amplification Using Polymerase Chain Reaction (PCR)

Name:

Group:

2. Mix the components briefly. Spin down the mixture by pulsing in a microcentrifuge to bring all reaction components to the bottom of the tube.

- 3. Aliquot b  $\mu$ L into each 0.2 ml PCR tube.
- 4. Add a µL of template DNA to each reaction except for the negative control tube (DNA will be replaced by 1.5 µL of ddH<sub>2</sub>O for the negative control).
- 5. Lastly, add 0.2 µL of Taq polymerase into each of 0.2 ml PCR tube.

#### 6. Carry out amplification in the thermocycler with the following PCR parameter.

#### Table 2.2: PCR parameter

Step	Temperature (°C)	Time	No. cycles
Pre-Denaturation	95		1
Denaturation	95		
Annealing			30
Extension	72		
Final Extension	72		1
Soak	4	00	00

7. After amplification of the products, visualize PCR products using 2% agarose

#### gel electrophoresis.

#### NOTES

- Make sure you collect the amplified products directly after the process reaches • the soaking step.
- Never let the products soaked for more than half an hour since this action might . cause damage to the thermocycler.

Name:	Matric No:	Group:
Results		[5 marks]



Practical 2: DNA Amplification Using Polymerase Chain Reaction (PCR)

.